

DISTRIBUTED BY:



# Crutches

The Physician's Choice crutches features durable construction, comfortable and adjustable padded hand grips and contoured underarm pads. The convenient push-button design allows you to quickly adjust the height in 1" increments. Non-skid tips provide excellent traction resulting in a stable performance. These crutches are available in Distribution in adult and tall sizes.



## FEATURES/BENEFITS:

- Aluminum crutch provides stability and durability (3/4 inch Frame diameter)
- Double-extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-sized vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Not made with natural rubber latex
- Aluminum construction provides support for up to 300 lbs.
- Comfortable, underarm design



MFG#	DESCRIPTION	PACKAGING	BRAND
SLI-FM-L-AL	Adult Tall Crutches, 5.9' - 6.5'	10 pair/case	Physician's Choice
SLI-FM-M-AL	Adult Crutches, 5.1' - 5.9'	10 pair/case	Physician's Choice

Contact your HMG Representative today for pricing and availability.



833-464-2426 | [info@hgm-distribution.com](mailto:info@hgm-distribution.com) | [www.hgm-distribution.com](http://www.hgm-distribution.com)