

DISTRIBUTED BY:



Electrosurgery

Electrosurgical Pencil & Pad

ELECTROSURGICAL PENCIL:

- CUT/yellow button for cutting.
- COAG/blue button for coagulation.
- Available with push button switch controls.
- Switches have a high level of tactile feedback surgeons can hear/feel the switch, preventing accidental activation.
- Slim, lightweight, tapered ergonomic streamlined design.
- Slip proof and water resistant.
- Buttons designed for smooth activation.
- Hexagonal shape to protect against inadvertent twisting.
- Electrodes lock safely into place and are easily replaced.



Electrosurgical Pad:

- Foam backing
- Adult model
- Split pad, with cable
- For Valleylab generator



MFG#	DESCRIPTION	PACKAGING	BRAND
18502450	Electrosurgical Pencil Button Switch w/ Non Stick Tip	100/cs	Forsure Medical
18507807	Electrosurgical Pad Adult w/ Cables	100/cs	Forsure Medical

Contact your HMG Representative today for pricing and availability.



833-464-2426 | info@hmg-distribution.com | www.hmg-distribution.com