

Preconnected Foley Catheterization Tray



PRODUCT DESCRIPTION:

The Preconnected Foley Catheterization Tray from Bridon Corp contains everything you need to perform a safe and effective catheterization procedure. This sterile kit provides clinicians flexibility to choose their preferred drainage or leg bag. Gloves and povidone-iodine swab sticks included in the kit help maintain sterility and reduce the risk of infection.

FOLEY CATHETERIZATION TRAY KIT INCLUDES:

- 2000cc Urine bag with anti-reflux valve
- 16 in. (40.64 cm) latex red urethral catheter, 2-way, 5cc
- 10cc pre-filled sterile water syringe
- CSR wrap, waterproof and fenestrated drapes
- (2) Nitrile powder-free gloves (large)
- (3) Povidone-iodine swab sticks
- Lubricating jelly 5 gram
- 2 oz. specimen container with patient label
- Paper box
- Plastic tray



ITEM#	DESCRIPTION	PACKAGING	BRAND
1232	16 Fr Preconnected Foley Cath Tray	1 ea/10 per case	Bridon Corp

Contact your HMG Representative today for pricing and availability.

