

DISTRIBUTED BY:



Blunt Fill Needles

Reli® Blunt Fill Needles eliminate the use of a sharp when reconstituting medication and reduce the risk of a needle-stick injury, while maintaining compatibility with most syringes available in the market.

FEATURES/BENEFITS:

- Angled bevel reduces coring of vial stopper
- 18G ensures rapid filling and smooth flow rate suitable for viscous medications
- 1.5" Needle aids penetration of different septum barriers
- Eliminates the use of a sharp, reducing the risk of a needle-stick injury
- Compatible with most syringes available in the market
- Sterile, Single use, PVC free, not made with natural rubber latex
- NOTE: RELI® Blunt Fill Needles are not intended for skin injection



MFG#	DESCRIPTION	PACKAGING	BRAND
F-BFN18G151	Blunt Fill Needle – Filtered, 18G X 1.5"	100EA/BX, 100BX/CS	MYCO

Contact your HMG Representative today for pricing and availability.



833-464-2426 | info@hmg-distribution.com | www.hmg-distribution.com