

# Safety Blood Collection Sets

Rospital's blood collection sets are designed to increase safety and efficiency when drawing blood samples. The easy to use luer adapter provides a safe collection environment for both clinician and patient. The luer adapter is already added to the butterfly needle therefore, the holder can immediately be screwed on without having to assemble the adapter first. In many cases this can save valuable time.

### FEATURES/BENEFITS:

- Double wings make puncture more safe
- Soft and transparent PVC tube can observe the vein blood flow clearly
- Luer Slip Adapter
- Sterile | EO
- Non-Toxic & Latex Free
- Safety Needle
- Non-pyrogenic
- Single-use only
- Ready to use right out of the package, with no assembly required
- 100 each per box – 10 boxes per case



MFG#	DESCRIPTION	PACKAGING	BRAND
SFYBC21	Safety Blood Collection Set without holder, 21g x 3/4"	1,000 ea/case	Rospital
SFYBC23	Safety Blood Collection Set without holder, 23g x 3/4"	1,000 ea/case	Rospital
SFYBC25	Safety Blood Collection Set without holder, 25g x 3/4"	1,000 ea/case	Rospital

Contact your HMG representative today for pricing and availability.

