

Luer Lock & Slip Tip Syringes

We make your life easier by providing the highest quality clinical syringes, such as slip tip and luer lock syringes, at the best prices. Luer Lock syringes are ideal for general use and for administering medications in a variety of settings. Not only do these syringes provide a secure connection between the injection needle and syringe, but they help prevent the needle from being detached inadvertently. Slip Tip syringes allow a needle to be quickly and conveniently pushed straight on to the end of the tip. They also help reduce dead space when using a large capacity syringe.



FEATURES/BENEFITS:

- 1 ml. Graduations
- Luer Lock and Slip Tip
- Excellent plunger slide characteristics
- Transparent barrel with clear graduations to allow exact dosing of fluid
- Secure backstop to prevent accidental plunger withdrawal
- Double seal stopper to help prevent medication leakage
- Large, stable finger grips to ensure ease of use
- Bold, precise scale markings make precise measurement of medication fast and easy
- Sterile - Single use only
- Non-toxic, Non-pyrogenic
- Not made with natural rubber latex



MFG#	DESCRIPTION	PACKAGING	BRAND
90100000	1mL Luer Lock Syringe	100 ea/innerbox, 32 innerboxes/case = 3,200 ea/case (single sterile)	Forsure Medical Products
90300000	3mL Luer Lock Syringe	100 ea/innerbox, 30 innerboxes/case = 3,000 ea/case (single sterile)	Forsure Medical Products
90300001	3mL Slip Tip Syringe	100 ea/innerbox, 30 innerboxes/case = 3,000 ea/case (single sterile)	Forsure Medical Products
90500000	5mL Luer Lock Syringe	100 ea/innerbox, 24 innerboxes/case = 2,400 ea/case (single sterile)	Forsure Medical Products
91000000	10mL Luer Lock Syringe	100 ea/innerbox, 16 innerboxes/case = 1,600 ea/case (single sterile)	Forsure Medical Products
91000001	10mL Slip Tip Syringe	100 ea/innerbox, 16 innerboxes/case = 1,600 ea/case (single sterile)	Forsure Medical Products
92000000	20 mL Slip Tip Syringe	50 ea/innerbox, 16 innerboxes/case = 800 ea/case (single sterile)	Forsure Medical Products
92000001	20cc Slip Tip Syringe	50 ea/innerbox, 16 innerboxes/case = 800 ea/case (single sterile)	Forsure Medical Products
93000000	30mL Luer Lock Syringe	50 ea/innerbox, 12 innerboxes/case = 600 ea/case (single sterile)	Forsure Medical Products
96000000	60mL Luer Lock Syringe	25 ea/innerbox, 16 innerboxes/case = 400 ea/case (single sterile)	Forsure Medical Products
96000003	60cc Catheter Tip Syringe	12 ea/innerbox, 25 innerboxes/case = 300 ea/case (single sterile)	Forsure Medical Products

Contact your HMG representative today for pricing and availability.

