

Co-Ease® COHESIVE



Co-Ease® cohesive bandages are designed to provide a consistent level of compression, that won't constrict or compromise proper circulation during extended wear. The permeable fabric is soft, comfortable and gentle on skin and features effortless bilateral tear ability. The self-adhering function allows the bandage to stick only to itself without the need for clips or fasteners. This lightweight bandage also provides strong, slip-resistant support and is water-resistant.



FEATURES/BENEFITS:

- Lightweight and self-adhesive bandage
- Extremely comfortable, non slip construction
- Made of soft, permeable fabric
- Water-resistant
- Effortless tear ability
- Provides strong and controlled support
- Offers long-term compression that won't constrict
- Available in multiple sizes

APPLICATIONS:

- Compression helps stop bleeding in trauma setting
- Provides long-term compression after blood donation
- Ideal for IV therapy and blood donor services
- Used as a fixation bandage for securing cold packs, braces, & padding
- Protects primary dressings
- Supplies a high-level of support to joint and limbs
- Conforms to the most difficult body contours
- Provides gentle wound care for geriatric patients

MFG#	DESCRIPTION	PACKAGING	BRAND
CE10LF	Cohesive Bandage, 1" x 5 yds, tan non latex	60 EA/BX & 12 BX/CS	Tempo
CE20LF	Cohesive Bandage, 2" x 5 yds, tan non latex	36 EA/BX & 12 BX/CS	Tempo
CE30LF	Cohesive Bandage, 3" x 5 yds, tan non latex	24 EA/BX & 12 BX/CS	Tempo
CE40LF	Cohesive Bandage, 4" x 5 yds, tan non latex	18 EA/BX & 12 BX/CS	Tempo
CE60LF	Cohesive Bandage, 6" x 5 yds, tan non latex	12 EA/BX & 12 BX/CS	Tempo

Contact your HMG representative today for pricing and availability.

